

Avalo® Medical Cart

Compatible with  **Quantra®**
Hemostasis Analyzer

Take your Quantra for a Ride!




- Rugged
- Stable
- Maneuverable
- Ample Storage
- Infection Prevention
- Safe and Secure

Avalo Medical Cart

Key Features and Benefits

- Configured for use with Quantra Hemostasis Analyzer
- Precision-welded steel frame with durable, high-impact panels
- Compatible with hospital-grade disinfectants
- Stable and secure with low center of gravity
- Security Bracket for Quantra
- Slim, lightweight design simplifies movement and steering
- Seamless drawer trays are easy to clean and open/close
- Wide selection of accessories and organization systems
- Ample storage accommodates:
 - Universal Power Supply
 - Thermal printer
 - Up to two cartridge kits
 - BD vacutainers
 - No-sharps container
 - Cleaning supplies
- Delivered fully assembled; no tools required
- 10-year limited warranty

Dimensions



Product Name	Part Number
Avalo Medical Cart with "HSQ" HemoSonics Configuration	AM8MC-LCD-N-DR011-HSQ

Request More Information

Avalo Medical Cart

Capsa Healthcare

jgrabowski@capsahealthcare.com

8170 Dove Parkway, Canal Winchester, OH 43110

708.274.0554

CapsaHealthcare.com

Quantra System

HemoSonics

contact@hemosonics.com

US: +1 (800) 280 5589

EU: +33 1 46 88 20 00

HemoSonics.com